



Earmuffs

SHORT DESCRIPTION

The hearing protection earmuffs are equipped with padded foam ear cushions and can be worn without any discomfort for long periods of time. The adjustable headband provides an unmatched comfort.

HOW DOES THIS TOOL WORK?

The earmuffs allow a noise reduction of 22 decibels, which facilitates focus on a task in a noisy context. In addition, children with auditory hyper-reactivity will better tolerate noisy environments. It is recommended to use the earmuffs with partial noise reduction, because very often it is not necessary to completely cut out the noise, but just alleviate it.

HOW TO USE IN EVERYDAY LIFE?

- They reduce unwanted parasitic noises (ex: other people talking, movements, street sounds).
- They foster and maintain attention during exams, individual work or during homework.
- They allow children with auditory hyper-reactivity (intolerance to sound and noises) to tolerate noisy activities (movie theater, concert, show, fireworks, vacuum cleaner, etc.)

WHO IS IT FOR?

The earmuffs are recommended for children, from 3 years old, who do not have a hearing problem (ex: tinnitus, auditory processing disorder, partial hearing loss). They can be suggested for children having difficulty functioning in a noisy environment because of sensory difficulties (hyperreactivity to noise) or attention difficulties. With the earmuffs children can also improve their focus during important tasks (reading, writing, exams, homework, etc.)

PRECAUTIONS

It is preferable to use the earmuffs following recommendation from a professional and for specific activities (reading, exams). A supervising adult could determine if the moment is appropriate to use the earmuffs, as some children will tend to use them during explanations or team work.